

CHO-HER2 | 305413MH**Description**

CHO-HER2

CHO

HER2

85,000

HER2

Organism**Tissue****Disease**

Chinese hamster ovary, non-neoplastic; genetically engineered for HER2 (ErbB2/CD340) surface expression (medium-high expression level)

Applications

Antibody screening; ADCC/CDC assays; HER2-targeted therapy development; breast/gastric cancer research; flow cytometry

Synonyms

CHO-HER2

Age**Gender****Morphology****Cell type**

Epithelial cells

Growth properties

/

Citation

CHO-HER2 High Cytion 305413H

Biosafety level

1

NCBI_TaxID

10029

Product sheet

CHO-HER2 | 305413MH

CellosaurusAccession CVCL_A8W7

GMO Status GMO-S1: This CHO derivative contains a medium-to-high HER2 expression construct for evaluating HER2-targeted therapeutics. This classification applies only within Germany and may differ elsewhere.

Receptors expressed HER2

Culture Medium A DMEM:Ham's F12 (1:1) w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO3 Cytion
InSCREENeX InSCREENeX INS-ME-1039

Supplements 5% FBS G418- 0.5 mg/mL

Dissociation Reagent -EDTA

Doubling time approx. 14-16 hours

Subculturing PBS PBS /EDTA T25 1 T75 3 37°C

Split ratio 1 to 5

Seeding density 2 to 5 x 10⁴ cells/cm²

Fluid renewal 2 3

Post-Thaw Recovery 1:2 1:3 T25 24

Freeze medium FBS +10% DMSO CM-1 Cytion 800100

Product sheet

CHO-HER2 | 305413MH

Thawing and
Culturing Cells

1.					
2.		-150°C			3
3.		37°C		40-60	
4.			70%		
5.			8	15	
6.	300 x g	3			
7.	10			T25	T25
8.					

Incubation
Atmosphere

37°C, 5% CO₂, humidified atmosphere.

Shipping
Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage
Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

Sterility

PCR