

Product sheet

CHO-PDCD1 | 305973

Description	/ 6,250		
	CHO-PDCD1	CHO	1 PD-1 PDCD1/CD279 PD-1
	CHO-PDCD1		-

Organism

Tissue			
--------	--	--	--

Disease	PD-1 PDCD1/CD279		
---------	------------------	--	--

Applications	PD-1	T	
--------------	------	---	--

Age			
-----	--	--	--

Gender

Morphology			
------------	--	--	--

Cell type

Growth properties	/		
-------------------	---	--	--

Citation	CHO-PDCD1	Cytion	305973
----------	-----------	--------	--------

Biosafety level	1		
-----------------	---	--	--

NCBI_TaxID	10029		
------------	-------	--	--

CellosaurusAccession	CVCL_A8X0		
----------------------	-----------	--	--

CHO-PDCD1 | 305973

GMO Status	GMO-S1	CHO	PDCD1									
Receptors expressed	PDCD1/CD279											
Culture Medium	DMEM:Ham's F12 (1:1) w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO3 Cytidine											
	CHO	A	InSCREENeX	InSCREENeX	INS-ME-1039							
Supplements	5% FBS		G418-	0.5 mg/mL								
Dissociation Reagent	-EDTA											
Doubling time	14-16											
Subculturing			PBS	PBS	5%	2-3	/EDTA	T25	1	T75	3	37
Split ratio	1 5											
Seeding density	2 5×10 ⁴ /											
Fluid renewal	2 3											
Post-Thaw Recovery	1:2	1:3	T25		24							
Freeze medium	FBS +10% DMSO		CM-1	Cytion	800100		CM-1		Cytion	800100		

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere

37°C, 5% CO₂,

Shipping Conditions

-78 °C

Storage Conditions

-150 -196 -80 °C

Sterility

PCR