

Product sheet

CHO-uPAR | 305978

Description	/ 6,250				
	CHO-uPAR	CHO	uPAR	PLAUR/CD87	uPAR
	CHO-uPAR		uPAR		GPI
					uPA-uPA

Organism

Tissue

Disease uPAR PLAUR/CD87

Applications	uPAR	/
--------------	------	---

Age

Gender

Morphology

Cell type

Growth properties	/
-------------------	---

Citation	CHO-UPAR Cytion	305978
----------	-----------------	--------

Biosafety level 1

NCBI_TaxID	10029
------------	-------

CellosaurusAccession CVCL_A8X4

CHO-uPAR | 305978

GMO Status	GMO-S1	CHO	PLAUR/uPAR
-------------------	--------	-----	------------

Surface antigens	uPAR	PLAUR/CD87
-------------------------	------	------------

Receptors expressed	TACD2	TROP2	GA733-1
----------------------------	-------	-------	---------

Culture Medium	DMEM:Ham's F12 (1:1)	w: 3.1 g/L	w: 2.5 mM L-	w: 15 mM HEPES	w: 0.5 mM	w: 1.2 g/L NaHCO ₃	Cytion
	CHO	A	InSCREENeX	InSCREENeX	INS-ME-1039		

Supplements	5% FBS	G418-	0.5 mg/mL
--------------------	--------	-------	-----------

Dissociation Reagent	-EDTA
-----------------------------	-------

Doubling time	14-16
----------------------	-------

Subculturing	PBS	PBS	5%	2-3	/EDTA	T25	1	T75	3	37
---------------------	-----	-----	----	-----	-------	-----	---	-----	---	----

Split ratio	1	5
--------------------	---	---

Seeding density	2	5×10 ⁴	/
------------------------	---	-------------------	---

Fluid renewal	2	3
----------------------	---	---

Post-Thaw Recovery	1:2	1:3	T25	24
---------------------------	-----	-----	-----	----

Freeze medium	FBS	+10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	-----	-----------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere

37°C, 5% CO₂,

Shipping Conditions

-78 °C

Storage Conditions

-150 -196 -80 °C

Sterility

PCR