

Product sheet

CHO-CD36 | 305979

Description	/ 6,250									
	CHO-CD36		CHO	CD36	CD36	B		IV GPIV	FAT	CD36
	CHO-CD36	-		CD36						

Organism

Tissue										
--------	--	--	--	--	--	--	--	--	--	--

Disease	CD36									
---------	------	--	--	--	--	--	--	--	--	--

Applications	CD36									
--------------	------	--	--	--	--	--	--	--	--	--

Age										
-----	--	--	--	--	--	--	--	--	--	--

Gender

Morphology										
------------	--	--	--	--	--	--	--	--	--	--

Cell type

Growth properties										
-------------------	--	--	--	--	--	--	--	--	--	--

Citation	CHO-CD36 Cytion 305979									
----------	------------------------	--	--	--	--	--	--	--	--	--

Biosafety level	1									
-----------------	---	--	--	--	--	--	--	--	--	--

NCBI_TaxID	10029									
------------	-------	--	--	--	--	--	--	--	--	--

CellosaurusAccession	CVCL_8848									
----------------------	-----------	--	--	--	--	--	--	--	--	--

Product sheet

CHO-CD36 | 305979

GMO Status	GMO-S1	CHO	CD36						
Receptors expressed	CD36								
Culture Medium	DMEM:Ham's F12 (1:1)	w: 3.1 g/L	w: 2.5 mM L-	w: 15 mM HEPES	w: 0.5 mM	w: 1.2 g/L NaHCO3			
	CHO	A	InSCREENeX	InSCREENeX	INS-ME-1039				
Supplements	5% FBS	G418-	0.5 mg/mL						
Dissociation Reagent	-EDTA								
Doubling time	14-16								
Subculturing	PBS	PBS	5%	2-3	/EDTA	T25	1	T75	3
Split ratio	1	5							
Seeding density	2	5×10 ⁴	/						
Fluid renewal	2	3							
Post-Thaw Recovery	1:2	1:3	T25	24					
Freeze medium	FBS	+10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100	

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere

37°C, 5% CO₂,

Shipping Conditions

-78 °C

Storage Conditions

-150 -196 -80 °C

Sterility

PCR