

Product sheet

HEK293-CD20 | 305987

Description	/ 6,250				
	HEK293-CD20	293 HEK293	CD20 MS4A1	B	CD20 B
	HEK293-CD20	CD20			

Organism

Tissue					
--------	--	--	--	--	--

Disease	/	HEK293			
---------	---	--------	--	--	--

Applications	CD20	CAR-T	ADCC/CDC	B	
--------------	------	-------	----------	---	--

Age					
-----	--	--	--	--	--

Gender

Morphology					
------------	--	--	--	--	--

Cell type

Growth properties					
-------------------	--	--	--	--	--

Citation	HEK293-CD20 Cytion	305987			
----------	--------------------	--------	--	--	--

Biosafety level	1				
-----------------	---	--	--	--	--

NCBI_TaxID	9606				
------------	------	--	--	--	--

CellosaurusAccession	CVCL_A8V4				
----------------------	-----------	--	--	--	--

HEK293-CD20 | 305987

GMO Status GMO-S1 HEK293 CD20 MS4A1

Receptors expressed CD20

Culture Medium RPMI 1640 w 2.0 mM w 2.0 g/L NaHCO₃ Cytion 820700a

Supplements 10% FBS 1 mM 10 mM HEPES 1% NEAA G418- 1 mg/mL

Dissociation Reagent -EDTA

Doubling time 24-36

Subculturing PBS PBS5% 2-3 /EDTA T25 1 T75 3 37°C

Split ratio 1 5

Seeding density 2 4×10^4 /

Fluid renewal 2 3

Post-Thaw Recovery 1:2 1:3 T25 24

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Product sheet

HEK293-CD20 | 305987

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5%_{CO2}

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR