

Product sheet

Hs 611.T | 305915

Description

Hs 611.T 48 cHL T cHL Hs 611.T  
Hs 611.T HRS NF-κB JAK-STAT CD30 CD15 CD25

Organism

Tissue

Disease

Metastatic site

Applications NF-κB JAK-STAT CD30/CD15/CD25

Synonyms Hs-611-T Hs611T HS611T

Age 48

Gender

Ethnicity

Morphology

Cell type B

Growth properties

Citation Hs 611.T Cytion 305915

Biosafety level 1

NCBI\_TaxID 9606

**Hs 611.T | 305915**

**CellosaurusAccession** CVCL\_0823

**Culture Medium** DMEM w 4.5 / w 4 L- w 3.7 / NaHCO<sub>3</sub> w 1.0 Cytion 820300a

**Supplements** 10% FBS

**Dissociation Reagent**

**Doubling time** 24 48

**Subculturing** 3–5 × 10<sup>5</sup> / 300×g 5

**Split ratio** 1 3

**Seeding density** 2 4 × 10<sup>5</sup> /

**Fluid renewal** 2 3

**Freeze medium** +10% DMSO

Product sheet

Hs 611.T | 305915

Thawing and  
Culturing Cells

- |    |         |        |     |       |
|----|---------|--------|-----|-------|
| 1. |         |        |     |       |
| 2. |         | -150°C |     | 3     |
| 3. |         | 37°C   |     | 40-60 |
| 4. |         |        | 70% |       |
| 5. |         | 8      | 15  |       |
| 6. | 200 x g | 5      |     |       |
| 7. |         |        |     |       |

Incubation  
Atmosphere

37°C, 5%<sub>CO2</sub>

Shipping  
Conditions

-78 °C

Storage  
Conditions

-150 -196 -80 °C