

P388-D1 | 400308

Description [P388 D1(IL-1)] -1 IL-1

Organism

Tissue

Disease

Synonyms P-388D1 P388D1 P388.D1 P3 88 D1

Breed/Subspecies DBA/2

Gender

Morphology

Cell type

Growth properties

Citation P388-D1 Cytion 400308

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_0477

Antigen expression H-2d

Product sheet

P388-D1 | 400308

Tumorigenic

Viruses	MAP B.piliformis	Ektromelie	K-	Kilham	Reo 3	PVM	LCM	M.pulmonis	MVM	Theiler's	GD VII	Toolan's	H-1	MH
---------	---------------------	------------	----	--------	-------	-----	-----	------------	-----	-----------	--------	----------	-----	----

Reverse transcriptase

MSI-status

Culture Medium	RPMI 1640	w 2.0 mM	w 2.0 g/L NaHCO3	Cytion	820700a
----------------	-----------	----------	------------------	--------	---------

Supplements	10% FBS
-------------	---------

Subculturing	5^5	/	5^5	1×10^6	/	3×10	5×10
--------------	-------	---	-------	-----------------	---	---------------	---------------

Seeding density	/	1×10^6
-----------------	---	-----------------

Fluid renewal

Post-Thaw Recovery	6	24	10
--------------------	-----	----	----

Freeze medium	FBS	+10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
---------------	-----	-----------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells

-
- 150°C
- 37°C
- 70%
- 8 15
- 300 x g 3
- 10 T25 T25
-

Incubation Atmosphere

37°C, 5% CO₂,

Flask Coating

Freezing Procedure

-78 °C

Shipping Conditions

-78 °C

Storage Conditions

-150 -196 -80 °C

Sterility

PCR