

SU-DHL-4 | 305106

[illegible]

SU-DHL-4 | 305106

Protein expression	IgG+ Kappa+ IgM- IgA- IgD- Lambda-	Bax Bak AIF	Caspase-9
---------------------------	------------------------------------	-------------	-----------

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO ₃ Cytion	820700a
-----------------------	--------------------	-------------------------------------	---------

Supplements	10% FBS
--------------------	---------

Doubling time	40
----------------------	----

Subculturing	⁵ / ⁵ 1×10 ⁶ / 3×10	5×10
---------------------	--	------

Split ratio	1:2 1:6
--------------------	---------

Fluid renewal	2 3
----------------------	-----

Freeze medium	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

Product sheet

SU-DHL-4 | 305106

Flask Coating

Freezing Procedure	-78 °C
--------------------	--------

Shipping Conditions	-78 °C
---------------------	--------

Storage Conditions	-150 -196 -80 °C
--------------------	----------------------

Sterility	PCR
-----------	-----