

HEK293T | 300189

[illegible]

HEK293T | 300189

Receptors expressed	Vitronectin													
Protein expression	CEA p53													
Tumorigenic														
Culture Medium	EMEM MEM Eagle w 2 mM L- w 2.2 g/L NaHCO3 w EBSS Cytion 820100a													
Supplements	10% FBS 1% NEAA													
Dissociation Reagent	Accutase													
Doubling time	30													
Subculturing		PBS	T25	3-5	PBS	T75	5-10	Accutase	T25	1-2	T75	2.5		
Seeding density	$1 \times 10^4 / \text{cm}^2$ 4													
Fluid renewal	2													
Post-Thaw Recovery														
Freeze medium		FBS	+10% DMSO	CM-1	Cytion		800100			CM-1	Cytion		800100	

Product sheet

HEK293T | 300189

Thawing and
Culturing Cells

1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
8.

Incubation
Atmosphere

37°C, 5% CO₂

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR