

P3X63Ag8 | 305171

Description	P3K27	BALB/c	MOPC-21	8-	20	/	HAT	3-
--------------------	-------	--------	---------	----	----	---	-----	----

Organism

Tissue	B
---------------	---

Synonyms	P3x63Ag8, P3-x63-Ag8, P3/x63-Ag8, P3/x63 Ag8, P3/x63/Ag8, P3-x63Ag8, P3x63 Ag8, P3x63 Ag8, P3 x 63Ag8, P3 x 63 Ag8, x63-Ag8, x63-AG8, x63-Ag8, P3x63, x63, GM03571
-----------------	--

Breed/Subspecies	BALB/c
-------------------------	--------

Gender**Morphology****Growth
properties**

Citation	P3x63Ag8	Cytion	305171
-----------------	----------	--------	--------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_3411
-----------------------------	-----------

**Protein
expression**

Antigen expression	H-2d
-------------------------------	------

Product sheet

P3X63Ag8 | 305171

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3	Cytion	820700a
----------------	--------------------	------------------	--------	---------

Supplements 10% FBS

Dissociation Reagent	Accutase
----------------------	----------

Doubling time 16 26

Subculturing	⁵ /	1×10
--------------	----------------	------

Split ratio 2×10⁵ /

Seeding density	3 5 × 10 ⁴ /m ²
-----------------	---------------------------------------

Fluid renewal 2 3

Freeze medium	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
---------------	---------------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Product sheet

P3X63Ag8 | 305171

**Incubation
Atmosphere** 37°C, 5% CO₂

Flask Coating

**Shipping
Conditions** -78 °C

**Storage
Conditions** -150 -196 -80 °C

Sterility PCR