

Y3-Ag 1.2.3 | 305207

Description	S210	8-	HAT	-	B	C. Milstein
--------------------	------	----	-----	---	---	-------------

Organism

Tissue

Disease

Synonyms	Y3.AG.1.2.3, Y3-Ag1.2.3, Y3-Ag1,2,3, Y3Ag1.2.3, Y-3-Ag 1.2.3, 210-RC Y3-Ag 1,2,3, 210RCY3-Ag1.2.3, 210RCY3-Ag123, Y3-Ag123, Y3, Y3M
-----------------	---

Breed/Subspecies	LOU/C
-------------------------	-------

Morphology

Growth properties

Citation	Y3-Ag 1.2.3 Cytion 305207
-----------------	---------------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	10116
-------------------	-------

CellosaurusAccession	CVCL_4342
-----------------------------	-----------

Protein expression

Y3-Ag 1.2.3 | 305207

Culture Medium DMEM w 4.5 / w 4 L- w 3.7 / NaHCO3 w 1.0 Cytion 820300a

Supplements 10% FBS

Subculturing ⁵ / 1×10

Fluid renewal 2 3

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

- Thawing and Culturing Cells**
- 1.
 2. -150°C 3
 3. 37°C 40-60
 4. 70%
 5. 8 15
 6. 300 x g 3
 7. 10 T25 T25
 - 8.

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Freezing Procedure -78 °C

Shipping Conditions -78 °C

Y3-Ag 1.2.3 | 305207

**Storage
Conditions**

-150 -196 -80 °C

/ /HLA

Sterility

PCR