

HK-CRISPR-CAP-D3-mEGFP | 301573

Description	HK-CRISPR-CAP-D3-mEGFP	Hela Kyoto	CRISPR/Cas9	mEGFP
		CAP-D3		

Organism

Tissue

Disease

Metastatic site

Applications	Condensin II	CAP-D3	CRISPR
---------------------	--------------	--------	--------

Synonyms	HK-CRISPR-CAP-D3-mEGFP #16, HK CRISPR CAP-D3-mEGFP
-----------------	--

Age	30
------------	----

Gender

Ethnicity

Morphology

Cell type

Growth properties

Citation	HK-CRISPR-CAP-D3-mEGFP Cytion 301573
-----------------	--------------------------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

HK-CRISPR-CAP-D3-mEGFP | 301573

CellosaurusAccession CVCL_UR44

Depositor EMBL

GMO Status GMO-S1 HeLa Kyoto CAP-D3 CRISPR mEGFP

Products EGFP

Culture Medium DMEM w 4.5 / w 4 L- w 3.7 / NaHCO3 w 1.0 Cytion 820300a

Supplements 10% FBS

Dissociation Reagent Accutase

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Fluid renewal 2 3

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

HK-CRISPR-CAP-D3-mEGFP | 301573

Thawing and Culturing Cells				
1.				
2.		-150°C		3
3.		37°C	40-60	
4.			70%	
5.		8	15	
6.	300 x g	3		
7.	10		T25	T25
8.				

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Freezing Procedure -78 °C

Shipping Conditions -78 °C

Storage Conditions -150 -196 -80 °C

/ /HLA

Sterility PCR