

Product sheet

HROC300 T2 M1 | 300866

Description	HROC300 T2 M1		HROC		"T2"		"M1"		
	HROC300 T2 M1								
	HROC				CIN	CpG	CIMP	MSS	MSI-H
Organism									
Tissue									
Disease	TNM	T4aN1bM1R2L0V1	G2 Lk(n) + 3 $\sum$ Lk(n) 22						
Metastatic site	TNM M1	T4aN1bM1R2							
Applications	MSS	HROC Linnebacher							
Age	73								
Gender									
Ethnicity									
Morphology									
Cell type									
Growth properties									
Citation	HROC300 T2 M1		Cytion	300866					
Biosafety level	1								
NCBI_TaxID	9606								

Product sheet

HROC300 T2 M1 | 300866

CellosaurusAccession CVCL_VQ94

GMO Status HROC Linnebacher

MSI-status MSS

Culture Medium DMEM:Ham's F12 1:1 w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO3 Cytion 8

Supplements 10% FBS

Dissociation Reagent Accutase

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Fluid renewal 3 5

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Product sheet

HROC300 T2 M1 | 300866

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 ° -80 °C

Sterility

PCR