

Product sheet

L Wnt-3A | 305184

Description	L Wnt-3A	L			Wnt	Wnt-3A	Wnt	
		L Wnt-3A	Wnt-3A		Wnt			Wnt
	Wnt-3A	L Wnt-3A		Wnt		Wnt		
Organism								
Tissue								
Synonyms	L-Wnt-3A, L-Wnt3A, LWnt3A, LWnt-3A							
Breed/Subspecies	C3H/An							
Age	100							
Gender								
Morphology								
Growth properties								
Citation	L Wnt-3A	Cytion	305184					
Biosafety level	1							
NCBI_TaxID	10090							
CellosaurusAccession	CVCL_0635							
GMO Status	GMO-S1	L	L Wnt-3A	PGK	Wnt3a		Wnt3a	L

L Wnt-3A | 305184

Protein expression Wnt-3A

Culture Medium DMEM w 4.5 / w 4 L- w 3.7 / NaHCO₃ w 1.0 Cytion 820300a

Supplements 10% FBS 0.4 mg/mL G-418

Dissociation Reagent Accutase

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Fluid renewal 2 3

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere 37°C, 5% CO₂

Product sheet

L Wnt-3A | 305184

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR