

Product sheet

HEK293T/17 | 305117

Description

293T/17 HEK293 SV40 T
293T/17 293T/17
293T/17

Organism

Tissue

Applications

Synonyms

HEK293T/17 HEK-293T/17 HEK 293T/17

Age

Gender

Morphology

Growth properties

Citation

HEK293T/17 Cytion 305117

Biosafety level

1

NCBI_TaxID

9606

CellosaurusAccession

CVCL_1926

GMO Status

GMO-S1 HEK293T/17 SV40

HEK293T/17 | 305117

Antigen expression

SV40 T

Viruses

SV40 SV40 T

Culture Medium
DMEM w 4.5 / w 4 L- w 3.7 / NaHCO₃ w 1.0 Cytion 820300a**Supplements**

10% FBS

Dissociation Reagent

Accutase

Subculturing

PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Fluid renewal

2 3

Freeze medium

FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Thawing and Culturing Cells

1.

2.

-150°C

3

3.

37°C

40-60

4.

70%

5.

8

15

6.

300 x g 3

7. 10

T25

T25

8.

Product sheet

HEK293T/17 | 305117

Incubation
Atmosphere 37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure -78 °C

Shipping
Conditions -78 °C

Storage
Conditions -150 -196 -80 °C

Sterility
PCR