

Product sheet

GH3 | 300383

Description

GH3
GH3
GH3
GH1
GH3
GH3

Organism

Tissue

Disease

Synonyms

GH 3

Breed/Subspecies

Wistar Furth

Age

7

Gender

Morphology

Growth properties

Citation

GH3 Cytion 300383

Biosafety level

1

NCBI_TaxID

10116

CellosaurusAccession

CVCL_0273

Products

Culture Medium	Ham's F12K	w 2.0 mM L-	w 2.0 mM	w 2.5 g/L NaHCO3	Cytion	820608a
-----------------------	------------	-------------	----------	------------------	--------	---------

Supplements	15%	2.5%	FBS
--------------------	-----	------	-----

Dissociation Reagent Accutase

Subculturing	15	PBS	T25	3-5	T75	5-10		Accutase	T25	1-2	T75	2.5
---------------------	----	-----	-----	-----	-----	------	--	----------	-----	-----	-----	-----

Freeze medium	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
---------------	---------------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere 37°C, 5% CO₂,

Flask Coating

Product sheet

GH3 | 300383

**Freezing
Procedure**

-78 °C

**Shipping
Conditions**

-78 °C

**Storage
Conditions**

-150 -196 -80 °C

Sterility

PCR