

Product sheet

CC531 | 500387

Description

CC531 --1,2- DMH Wistar CC531
CC531

Organism

Tissue

Disease

Synonyms

CC-531

Breed/Subspecies

WAG

Growth
properties

Citation

CC531 Cytion 500387

Biosafety level

1

NCBI_TaxID

10116

CellosaurusAccession

CVCL_0206

Tumorigenic

WAG-Rij

Culture
Medium

RPMI 1640 w 2.0 mM w 2.0 g/L NaHCO3 Cytion 820700a

CC531 | 500387

Supplements 10% FBS 20 mM HEPES

Dissociation Reagent Accutase

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Seeding density $1 \times 10^4 / \text{cm}^2$ 3 4

Fluid renewal 2 3

Post-Thaw Recovery 4 5×10^4 48

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere 37°C, 5% CO_2

Flask Coating

Product sheet

CC531 | 500387

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR