

Product sheet

HROC147 T0 M1 | 300856

Description - 2006 CRC

Organism

Tissue UICC IV CRC TNM T3N2M1R0L1V1 G3 Lk(n) +4 Σ Lk(n) 32

Disease

Synonyms HROC147

Age 54

Gender

Ethnicity

Morphology

Growth properties

Citation HROC147 T0 M1 Cytion 300856

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_1G09

Protein expression PTEN

HROC147 T0 M1 | 300856

Tumorigenic**Viruses**

SV40 JC/BK HBV HCV HIV

MSI-status

MSI-L

Mutational profile

APCmut p53wt K-Rasmut B-RAFwt N-Raswt H-Raswt PIK3CAwt

Culture MediumDMEM:Ham's F12 1:1 w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO₃ Cytion 8**Supplements**

10% FBS

Dissociation Reagent

Accutase

Doubling time

29

Subculturing

PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Seeding density 2×10^4 /**Fluid renewal**

3 5

Post-Thaw Recovery**Freeze medium**

FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Product sheet

HROC147 T0 M1 | 300856

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR