

EB3 | 300373

Description

EB3

EB3

Organism

Tissue

Disease

Metastatic site

Applications

Synonyms

EB-3 Epstein-Barr-3 GM04679

Age

3

Gender

Ethnicity

Morphology

Cell type

B

Growth properties

Citation

EB3 Cytion 300373

Biosafety level

2

NCBI_TaxID

9606

EB3 | 300373

CellosaurusAccession CVCL_1185

Surface antigens HLA A3 Aw32 Cw2

Isoenzymes G6PD, A

Viruses EB EBNA

Culture Medium RPMI 1640 w 2.0 mM w 2.0 g/L NaHCO₃ Cytion 820700a

Supplements 10% FBS

Subculturing ⁵ / 1×10

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

EB3 | 300373

**Incubation
Atmosphere** 37°C, 5% CO₂

Flask Coating

**Freezing
Procedure** -78 °C

**Shipping
Conditions** -78 °C

**Storage
Conditions** -150 -196 -80 °C

/ **/HLA**

Sterility PCR