

Product sheet

OS-RC-2 | 305086

|             |         |     |          |         |       |         |     |  |
|-------------|---------|-----|----------|---------|-------|---------|-----|--|
| Description | OS-RC-2 | RCC | RCC      | RCC     |       |         |     |  |
|             | OS-RC-2 |     |          | OS-RC-2 | 40    | 2       | 1   |  |
|             | OS-RC-2 | RCC | -α TNF-α | -6 IL-6 | TNF-α | OS-RC-2 | DNA |  |

Organism

Tissue

Disease

|          |             |
|----------|-------------|
| Synonyms | OSRC2, RC-2 |
|----------|-------------|

|     |    |
|-----|----|
| Age | 52 |
|-----|----|

Gender

|           |  |
|-----------|--|
| Ethnicity |  |
|-----------|--|

Morphology

|                   |  |
|-------------------|--|
| Growth properties |  |
|-------------------|--|

|          |                       |
|----------|-----------------------|
| Citation | OS-RC-2 Cytion 305086 |
|----------|-----------------------|

|                 |   |
|-----------------|---|
| Biosafety level | 1 |
|-----------------|---|

|            |      |
|------------|------|
| NCBI_TaxID | 9606 |
|------------|------|

|                      |           |
|----------------------|-----------|
| CellosaurusAccession | CVCL_1626 |
|----------------------|-----------|

Product sheet

OS-RC-2 | 305086

Tumorigenic

Culture Medium

RPMI 1640 w 2.0 mM w 2.0 g/L NaHCO3 Cytion 820700a

Supplements

10% FBS

Dissociation Reagent

Accutase

Subculturing

PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Split ratio

1:2 1:4

Fluid renewal

2 3

Freeze medium

FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere

37°C, 5% CO<sub>2</sub>

Product sheet

OS-RC-2 | 305086

Flask Coating

|                    |        |
|--------------------|--------|
| Freezing Procedure | -78 °C |
|--------------------|--------|

|                     |        |
|---------------------|--------|
| Shipping Conditions | -78 °C |
|---------------------|--------|

|                    |                      |
|--------------------|----------------------|
| Storage Conditions | -150   -196   -80 °C |
|--------------------|----------------------|

|           |     |
|-----------|-----|
| Sterility | PCR |
|-----------|-----|