

FAMPAC-A | 300310**Description** FAMPAC 43**Organism****Tissue****Disease****Synonyms** PA-CLS-13A PA-CLS 13A**Age** 43**Gender****Ethnicity****Morphology****Growth
properties****Citation** FAMPAC-A Cytion 300310**Biosafety level** 1**NCBI_TaxID** 9606**CellosaurusAccession** CVCL_5750**Antigen
expression** FAMPAC-A 12 Kras GGT(Gly)>GTT(Val)

Product sheet

FAMPAC-A | 300310

Tumorigenic

Karyotype

Culture Medium	DMEM:Ham's F12	1:1	w: 3.1 g/L	w: 2.5 mM L-	w: 15 mM HEPES	w: 0.5 mM	w: 1.2 g/L NaHCO3	Cytion	8
----------------	----------------	-----	------------	--------------	----------------	-----------	-------------------	--------	---

Supplements	10% FBS
-------------	---------

Dissociation Reagent	Accutase
----------------------	----------

Doubling time	32
---------------	----

Subculturing	PBS	T25	3-5	PBS	T75	5-10	Accutase	T25	1-2	T75	2.5
--------------	-----	-----	-----	-----	-----	------	----------	-----	-----	-----	-----

Fluid renewal	2	3
---------------	---	---

Post-Thaw Recovery	4	5×10	24
--------------------	---	------	----

Freeze medium	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
---------------	---------------	-------------	--------	-------------	--------

Product sheet

FAMPAC-A | 300310

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | | 8 | 15 | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR