

Product sheet

MKN-45 | 300489

Description

MKN-45

MKN-45

MKN-45

CEA E-cadherin

MKN-45

Organism

Tissue

Disease

Metastatic site

Synonyms

MKN 45, MKN45

Age

62

Gender

Ethnicity

Growth properties

/

Citation

MKN-45 Cytion 300489

Biosafety level

1

NCBI_TaxID

9606

CellosaurusAccession

CVCL_0434

MKN-45 | 300489

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	--------------------	------------------------------	--------	---------

Supplements	20%	FBS
--------------------	-----	-----

Dissociation Reagent	Accutase
-----------------------------	----------

Subculturing	15	PBS	T25	3-5	T75	5-10	Accutase	T25	1-2	T75	2.5
---------------------	----	-----	-----	-----	-----	------	----------	-----	-----	-----	-----

Freeze medium	FBS	+10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	-----	-----------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

Flask Coating

Freezing Procedure	-78 °C
---------------------------	--------

Product sheet

MKN-45 | 300489

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR