

E11 | 400494

Description	E11	E11	SV40 T	IFN	H-2kb
	E11	40	E11		
	E11		--		podocin CD2AP

Organism

Tissue

Breed/Subspecies	(CBA/Ca x C57BL/10)Tg(H2KbtsA58)
-------------------------	----------------------------------

Age

Gender

Cell type

Growth properties

Citation	E11 Cytion	400494
-----------------	------------	--------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_5737
-----------------------------	-----------

GMO Status	GMO-S1	Immorto	SV40 T
-------------------	--------	---------	--------

E11 | 400494

Protein expression	WT1 Lmx1b	NEPH1 FAT P-cadherin CD2AP ZO-1 podocalyxin podoplanin synpo podocin TRPC6 GAPDH
---------------------------	-----------	--

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3 Cytion	820700a
-----------------------	--------------------	-------------------------	---------

Supplements	10% FBS
--------------------	---------

Dissociation Reagent	Accutase
-----------------------------	----------

Subculturing	PBS	T25	3-5	PBS T75	5-10	Accutase	T25	1-2	T75	2.5
---------------------	-----	-----	-----	---------	------	----------	-----	-----	-----	-----

Seeding density	T75 ⁴	/	1×10 ³	33 /5% CO ₂	75%
------------------------	------------------	---	-------------------	------------------------	-----

Fluid renewal	2	3
----------------------	---	---

Freeze medium	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

Thawing and Culturing Cells	1.				
	2.		-150°C		3
	3.		37°C		40-60
	4.			70%	
	5.		8	15	
	6.	300 x g	3		
	7.	10		T25	T25
	8.				

Product sheet

E11 | 400494

**Incubation
Atmosphere** 33°C, 5% CO₂

Flask Coating

**Freezing
Procedure** -78 °C

**Shipping
Conditions** -78 °C

**Storage
Conditions** -150 -196 -80 °C

/ **/HLA**

Sterility PCR