

Product sheet

KHM-5M | 305148

Description

KHM-5M

8 IL-8

KHM-5M

KHM-5M

Organism

Tissue

Disease

Metastatic site

Synonyms

KHM/5M, KHM5M

Age

65

Gender

Morphology

Growth  
properties

Citation

KHM-5M Cytion 305148

Biosafety level

1

NCBI\_TaxID

9606

CellosaurusAccession

CVCL\_2975

Product sheet

KHM-5M | 305148

|                |                    |                  |        |         |
|----------------|--------------------|------------------|--------|---------|
| Culture Medium | RPMI 1640 w 2.0 mM | w 2.0 g/L NaHCO3 | Cytion | 820700a |
|----------------|--------------------|------------------|--------|---------|

Supplements 10% FBS

|                      |          |
|----------------------|----------|
| Dissociation Reagent | Accutase |
|----------------------|----------|

Doubling time 27

|              |     |     |     |     |     |      |          |     |     |     |     |
|--------------|-----|-----|-----|-----|-----|------|----------|-----|-----|-----|-----|
| Subculturing | PBS | T25 | 3-5 | PBS | T75 | 5-10 | Accutase | T25 | 1-2 | T75 | 2.5 |
|--------------|-----|-----|-----|-----|-----|------|----------|-----|-----|-----|-----|

Fluid renewal 2 3

|               |               |      |        |        |      |        |        |
|---------------|---------------|------|--------|--------|------|--------|--------|
| Freeze medium | FBS +10% DMSO | CM-1 | Cytion | 800100 | CM-1 | Cytion | 800100 |
|---------------|---------------|------|--------|--------|------|--------|--------|

Thawing and Culturing Cells

1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
8.

|                       |               |
|-----------------------|---------------|
| Incubation Atmosphere | 37°C, 5% CO2' |
|-----------------------|---------------|

Flask Coating

Product sheet

KHM-5M | 305148

|                    |        |
|--------------------|--------|
| Freezing Procedure | -78 °C |
|--------------------|--------|

|                     |        |
|---------------------|--------|
| Shipping Conditions | -78 °C |
|---------------------|--------|

|                    |                      |
|--------------------|----------------------|
| Storage Conditions | -150   -196   -80 °C |
|--------------------|----------------------|

|           |     |
|-----------|-----|
| Sterility | PCR |
|-----------|-----|