

SK-MEL-28 | 300337

Description T.Takahashi 51**Organism****Tissue****Disease****Synonyms** SK-Mel-28, SK.MEL.28, SK-MEL 28, SK MEL-28, SK MEL 28, SK Mel 28, SKMel-28, SKMEL-28, SK-MEL28, SK-Mel28, SK Mel28, SKMEL28, SKMel28, SKmel28, SKML-28, SK28, AU-Mel, P-36, P36**Age** 51**Gender****Morphology****Growth
properties****Citation** SK-MEL-28 Cytion 300337**Biosafety level** 1**NCBI_TaxID** 9606**CellosaurusAccession** CVCL_0526**Protein
expression** P53**Isoenzymes** PGM3 1 PGM1 1 ES-D 1 AK-1 1-2 GLO-1 2 G6PD B

SK-MEL-28 | 300337

Tumorigenic**Products****Mutational profile**

BRAF V600E DNA RT-PCR Western V600E BRAF

Culture MediumDMEM w 4.5 / w 4 L- w 3.7 / NaHCO₃ w 1.0 Cytion 820300a**Supplements**

10% FBS

Dissociation Reagent

Accutase

Subculturing

PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Split ratio

1:3

Fluid renewal

2 3

Post-Thaw Recovery

48

Freeze medium

FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Product sheet

SK-MEL-28 | 300337

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5% CO₂

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR

SK-MEL-28 | 300337

STR

Amelogenin: x,y
CSF1PO: 10,12
D13S317: 11,12
D16S539: 9,12
D5S818: 13
D7S820: 10
TH01: 7
TPOX: 8,12
vWA: 16,19
D3S1358: 16,18
D21S11: 28 29
D18S51: 12,16
Penta E: 8,12
Penta D: 9,1
D8S1179: 13
FGA: 19

HLA

A*: '11:01:01
B*: '40:01:02
C*: '03:04:01
DRB1*: '04:04:01
DQA1*: '03:01:01
DQB1*: '03:02:01
DPB1*: '03:01:01
E: '01:03:02