

Product sheet

M-MSV-Balb/3T3 | 400458

Description	M-MSV-Balb/3T3	BALB/c	3T3	3T3
	M-MSV		M-MSV-Balb/3T3	
Organism				
Tissue				
Synonyms	M-MSV-BALB/3T3			
Breed/Subspecies	BALB/c			
Age	14	17		
Gender				
Morphology				
Cell type				
Growth properties				
Citation	M-MSV-Balb/3T3	Cytion	400458	
Biosafety level	1			
NCBI_TaxID	10090			
CellosaurusAccession	CVCL_5793			
GMO Status	GMO-S1	M-MSV-Balb/3T3	MOMSV	Balb/3T3

M-MSV-Balb/3T3 | 400458

**Antigen
expression**

H-2d

Tumorigenic**Viruses**
**Reverse
transcriptase**
**Culture
Medium**
DMEM w 4.5 / w 4 L- w 3.7 / NaHCO₃ w 1.0

Cytion 820300a

Supplements

10% FBS

**Dissociation
Reagent**

Accutase

Subculturing

PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

**Seeding
density**
0.7 1×10^6 /**Fluid renewal**

2 3

**Freeze
medium**

FBS +10% DMSO CM-1 Cytion 800100

CM-1 Cytion 800100

Product sheet

M-MSV-Balb/3T3 | 400458

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5% CO₂

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR