

Product sheet

HK EGFP-LaminB1/H2B-mCherry | 300919

Description	HK EGFP-LaminB1/H2B-mCherry	Hela	EGFP
	mCherryH2B	DNA	DNAH2BmCherry
Organism			
Tissue			
Disease			
Metastatic site			
Applications	Lamin B1	H2B	
Synonyms	HeLa Kyoto EGFP-LaminB1H2B-mCherry		
Age	30		
Gender			
Ethnicity			
Morphology			
Cell type			
Growth properties			
Citation	HK EGFP-LaminB1/H2B-mCherry	Cytion	300919
Biosafety level	1		
NCBI_TaxID	9606		

HK EGFP-LaminB1/H2B-mCherry | 300919

CellosaurusAccession CVCL_UR41

Depositor EMBL

GMO Status GMO-S1 HeLa EGFP- B1 H2B-mCherry

Protein expression EGFP-LaminB1/H2B-mCherry

Products H2B

Culture Medium DMEM w 4.5 / w 4 L- w 3.7 / NaHCO₃ w 1.0 Cytion 820300a

Supplements 10% FBS

Dissociation Reagent Accutase

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Seeding density 1×10^4 /

Fluid renewal 2 3

Post-Thaw Recovery 4×10^4 5×10^4 24

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Product sheet

HK EGFP-LaminB1/H2B-mCherry | 300919

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | | 8 | 15 | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5%_{CO2}

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

Sterility

PCR