

Ba/F3 | 305224

Description	BA/F3	BALB/c	B	BaF3					
	Baf3	IL-3		Ba/F3	F3	Ba/F3	F3		IL-3
		BA/F3				CML	BCR-ABL	Ba/F3	
	BA/F3								

TKIs

Organism	
Tissue	

Synonyms	BA/F3, BaF3, BAF3, Baf3
-----------------	-------------------------

Breed/Subspecies	C3H
-------------------------	-----

Morphology	
Cell type	B

Growth properties	
--------------------------	--

Citation	Ba/F3 Cytion 305224
-----------------	---------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_0161
-----------------------------	-----------

Karyotype	Ba/F3	33%
------------------	-------	-----

Ba/F3 | 305224

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3	Cytion	820700a
-----------------------	--------------------	------------------	--------	---------

Supplements	5%	FBS	10 ng/mL	IL-3
--------------------	----	-----	----------	------

Subculturing	5^5	1×10^6	3×10^5	5×10^5
---------------------	-------	-----------------	-----------------	-----------------

Freeze medium	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	---------------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells	1.				
	2.		-150°C		3
	3.		37°C	40-60	
	4.		70%		
	5.		8	15	
	6.	300 x g	3		
	7.	10		T25	T25
	8.				

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

Flask Coating	
----------------------	--

Freezing Procedure	-78 °C
---------------------------	--------

Shipping Conditions	-78 °C
----------------------------	--------

