

DI TNC1 | 305343

GMO Status	GMO-S1	DI TNC1	GFAP	SV40
-------------------	--------	---------	------	------

Protein expression	$\alpha 2$
---------------------------	------------

Tumorigenic

Viruses	Simian virus 40 (SV40)
----------------	------------------------

Culture Medium	DMEM w 4.5 / w 4 L- w 3.7 / NaHCO3 w 1.0	Cytion 820300a
-----------------------	--	----------------

Supplements 10% FBS

Dissociation Reagent	Accutase
-----------------------------	----------

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Fluid renewal	2 3
----------------------	-----

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

DI TNC1 | 305343

Thawing and Culturing Cells				
1.				
2.		-150°C		3
3.		37°C	40-60	
4.			70%	
5.		8	15	
6.	300 x g	3		
7.	10		T25	T25
8.				

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Freezing Procedure -78 °C

Shipping Conditions -78 °C

Storage Conditions -150 -196 -80 °C

/ /HLA

Sterility PCR