

B-LCL-CDG3 | 302014

Description	B-LCL-CDG3	EBV	B	PMM2-CDG	PMM2*	CDG	PMM2	2 phosphomannomut
	B-LCL-CDG3	EBV-immortalized B		*PMM2*			PMM2	

Organism

Tissue

Disease

Applications	CDG	B	HLA
---------------------	-----	---	-----

Gender

Ethnicity

Morphology

Cell type	B
------------------	---

Growth properties

Citation	B-LCL-CDG3 Cytion 302014
-----------------	--------------------------

Biosafety level	2
------------------------	---

NCBI_TaxID	9606
-------------------	------

Viruses	EBV
----------------	-----

B-LCL-CDG3 | 302014

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	--------------------	------------------------------	--------	---------

Supplements	10%	FBS		
--------------------	-----	-----	--	--

Subculturing	⁵ /	2×10 ⁵ /	1×10 ⁵	5×10
---------------------	----------------	---------------------	-------------------	------

Fluid renewal

Post-Thaw Recovery				
---------------------------	--	--	--	--

Freeze medium	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

Thawing and Culturing Cells	1.			
	2.	-150°C		3
	3.	37°C	40-60	
	4.		70%	
	5.	8	15	
	6.	300 x g	3	
	7.	10	T25	T25
	8.			

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

Flask Coating	
----------------------	--

B-LCL-CDG3 | 302014

Freezing Procedure -78 °C

Shipping Conditions -78 °C

Storage Conditions -150 -196 -80 °C

/ **/HLA**

Sterility PCR