

Product sheet

CHO-CXCR4 | 305411L

Description

CHO-CXCR4-Medium-high CHO CXCR4 9500
CXCR4 CD184 HIV CXCL12

Organism

Tissue

Synonyms CHO-CXCR4

Age

Gender

Morphology

Growth properties /

Citation CHO-CXCR4 Medium-high Cytion 305411MH

Biosafety level 1

NCBI_TaxID 10029

Product sheet

CHO-CXCR4 | 305411L

Shipping
Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage
Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

/ /HLA

Sterility

PCR