

WM-115 | 305457

| | | | |
|--------------------|--------|--------|--------|
| Description | WM-115 | VGP | WM-115 |
| | WM-115 | | WM-115 |
| | WM-115 | WM-115 | |

Organism

| |
|---------------|
| Tissue |
|---------------|

Disease

| |
|------------------------|
| Metastatic site |
|------------------------|

Synonyms WM-115, WM 115, WM115F, WM115-mel, WM115mel, WC00079

| | |
|------------|----|
| Age | 55 |
|------------|----|

Gender

| |
|------------------|
| Ethnicity |
|------------------|

Growth properties

| | | |
|-----------------|-------|--------|
| Citation | WM115 | 305457 |
|-----------------|-------|--------|

Biosafety level 1

| | |
|-------------------|------|
| NCBI_TaxID | 9606 |
|-------------------|------|

CellosaurusAccession CVCL_0040

WM-115 | 305457

| | |
|---------------------------|-------------|
| Mutational profile | p.Val600Asp |
|---------------------------|-------------|

| | |
|-----------------------|-----------------------------------------------------------------|
| Culture Medium | EMEM MEM Eagle w 2 mM L- w 2.2 g/L NaHCO3 w EBSS Cytion 820100a |
|-----------------------|-----------------------------------------------------------------|

| | |
|--------------------|--------|
| Supplements | 10% 1% |
|--------------------|--------|

| | |
|-----------------------------|----------|
| Dissociation Reagent | Accutase |
|-----------------------------|----------|

| | |
|------------------------|-----------------------|
| Seeding density | 1 3x10 ⁴ / |
|------------------------|-----------------------|

| | |
|----------------------|-----------|
| Freeze medium | +10% DMSO |
|----------------------|-----------|

| | | | |
|------------------------------------|----|---------|-------|
| Thawing and Culturing Cells | 1. | | |
| | 2. | -150°C | 3 |
| | 3. | 37°C | 40-60 |
| | 4. | 70% | |
| | 5. | 8 | 15 |
| | 6. | 200 x g | 5 |
| | 7. | | |

| | |
|------------------------------|--------------------------|
| Incubation Atmosphere | 37°C, 5% CO ₂ |
|------------------------------|--------------------------|

| | |
|----------------------|--|
| Flask Coating | |
|----------------------|--|

WM-115 | 305457

**Shipping
Conditions**

-78 °C

**Storage
Conditions**

-150 -196 -80 °C

/ /HLA