



Sf9 | 604329

**Virus susceptibility** Baculoviruses, Autographa californica (MNPV), St.

**Culture Medium** PAN

**Supplements** 2% FBS

**Dissociation Reagent** Accutase

**Subculturing** 15 5 300xg 3

**Seeding density**  $1 \times 10^4$  cells/cm<sup>2</sup> 26-30

**Fluid renewal** 2 3

**Freeze medium** FBS +10% DMSO CM-1 Cytion 800100

Sf9 | 604329

**Thawing and Culturing Cells**

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

**Incubation Atmosphere**

27°C, 0% CO<sub>2</sub>, humidified atmosphere.

**Shipping Conditions**

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

**Storage Conditions**

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

/ /HLA

**Sterility**

PCR