



**NUGC-4 | 305645**

<b>Culture Medium</b>	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3 Cytion	820700a
-----------------------	--------------------	-------------------------	---------

<b>Supplements</b>	10% FBS
--------------------	---------

<b>Dissociation Reagent</b>	Accutase
-----------------------------	----------

<b>Doubling time</b>	29.9
----------------------	------

<b>Seeding density</b>	1 $4 \times 10^4 / \text{cm}^2$
------------------------	---------------------------------

<b>Fluid renewal</b>	2 3
----------------------	-----

<b>Freeze medium</b>	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

<b>Thawing and Culturing Cells</b>	1.			
	2.	-150°C		3
	3.	37°C		40-60
	4.		70%	
	5.	8	15	
	6.	300 x g	3	
	7.	10		T25 T25
	8.			

<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>
------------------------------	--------------------------

**NUGC-4 | 305645**

---

**Shipping  
Conditions**

-78 °C

**Storage  
Conditions**

-150 -196 -80 °C

/ **/HLA**

**Sterility**

PCR