

HEK293-PSMA | 305992

<b>Description</b>	/	6,250							
	HEK293-PSMA	293 HEK293	PSMA	II FOLH1/GCPII	PSMA				II
	HEK293-PSMA	PSMA	T	CAR-T					

**Organism**

**Tissue**

**Age**

**Gender**

**Morphology**

**Growth properties**

<b>Citation</b>	HEK293-PSMA Cytion 305992
-----------------	---------------------------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>Receptors expressed</b>	PSMA
----------------------------	------

**Product sheet**



**HEK293-PSMA | 305992**

<b>Culture Medium</b>	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO <sub>3</sub>	Cytion	820700a
<b>Supplements</b>	10% FBS 1 mM	10 mM HEPES	1% NEAA	G418- 1 mg/mL
<b>Dissociation Reagent</b>	-EDTA			
<b>Subculturing</b>		PBS	PBS5%	2-3 /EDTA T25 1 T75 3 37°C
<b>Fluid renewal</b>	2 3			
<b>Post-Thaw Recovery</b>	1:2 1:3	T25	24	
<b>Freeze medium</b>	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion 800100
<b>Thawing and Culturing Cells</b>	1.			
	2.	-150°C		3
	3.	37°C	40-60	
	4.		70%	
	5.	8	15	
	6.	300 x g	3	
	7.	10	T25	T25
	8.			
<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>			

HEK293-PSMA | 305992

---

**Shipping  
Conditions**

-78 °C

**Storage  
Conditions**

-150 -196 -80 °C

/ /HLA

**Sterility**

PCR