



**B-LCL-HROC32 (Bc HROC32) | 300862**

<b>Surface antigens</b>	CD19
-------------------------	------

<b>Viruses</b>	EBV
----------------	-----

<b>Culture Medium</b>	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3	Cytion	820700a
-----------------------	--------------------	------------------	--------	---------

<b>Supplements</b>	10%	FBS
--------------------	-----	-----

<b>Subculturing</b>	5	/	1x10
---------------------	---	---	------

<b>Freeze medium</b>	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	---------------	------	--------	--------	------	--------	--------

<b>Thawing and Culturing Cells</b>	1.				
	2.		-150°C	3	
	3.		37°C	40-60	
	4.		70%		
	5.		8	15	
	6.	300 x g	3		
	7.	10		T25	T25
	8.				

<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>
------------------------------	--------------------------

<b>Flask Coating</b>
----------------------

**B-LCL-HROC32 (Bc HROC32) | 300862**

**Freezing Procedure** -78 °C

**Shipping Conditions** -78 °C

**Storage Conditions** -150 -196 -80 °C

/ **/HLA**

**Sterility** PCR

**HLA**  
**A\***: '01:01:01, '02:01:01  
**B\***: '08:01:01, '51:01:01  
**C\***: '07:01:01, '14:02:01  
**DRB1\***: '03:01:01, '11:01:01  
**DQA1\***: '05:01:01, '05:05:01  
**DQB1\***: '02:01:01, '03:01:01  
**DPB1\***: '01:01:01, '04:01:01  
**E**: '01:01, '01:03