

B-LCL-HROC32 (Bc HROC32) | 300862

Surface antigens	CD19
-------------------------	------

Viruses	EBV
----------------	-----

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3	Cytion	820700a
-----------------------	--------------------	------------------	--------	---------

Supplements	10%	FBS
Subculturing	⁵ /	1×10

Freeze medium	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	---------------	------	--------	--------	------	--------	--------

Thawing and Culturing Cells	1.			
	2.	-150°C		3
	3.	37°C		40-60
	4.		70%	
	5.	8	15	
	6.	300 x g	3	
	7.	10		T25
	8.			T25

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

B-LCL-HROC32 (Bc HROC32) | 300862

Flask Coating

Freezing Procedure

-78 °C

Shipping Conditions

-78 °C

Storage Conditions

-150 -196 -80 °C

/ **/HLA**

Sterility

PCR

HLA

A*: '01:01:01, '02:01:01
B*: '08:01:01, '51:01:01
C*: '07:01:01, '14:02:01
DRB1*: '03:01:01, '11:01:01
DQA1*: '05:01:01, '05:05:01
DQB1*: '02:01:01, '03:01:01
DPB1*: '01:01:01, '04:01:01
E: '01:01, '01:03