



**B-LCL-HROC69 | 300864**

<b>Surface antigens</b>	CD19
-------------------------	------

<b>Viruses</b>	EBV
----------------	-----

<b>Culture Medium</b>	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO <sub>3</sub> Cytion	820700a
-----------------------	--------------------	-------------------------------------	---------

<b>Supplements</b>	10%	FBS
<b>Subculturing</b>	<sup>5</sup> /	1×10

<b>Freeze medium</b>	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

<b>Thawing and Culturing Cells</b>	1.			
	2.	-150°C		3
	3.	37°C		40-60
	4.		70%	
	5.	8		15
	6.	300 x g	3	
	7.	10		T25
	8.			T25

<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>
------------------------------	--------------------------

**B-LCL-HROC69 | 300864**

**Flask Coating**

**Freezing  
Procedure**

-78 °C

**Shipping  
Conditions**

-78 °C

**Storage  
Conditions**

-150 -196 -80 °C

/ **/HLA**

**Sterility**

PCR