

WEHI-3 | 400381

Description	WEHI-3	BALB/c
	WEHI-3	
Organism		
Tissue		
Disease		
Synonyms	WEHI 3, WEHI3, WEHi-3	
Breed/Subspecies	BALB/c	
Morphology		
Cell type		
Growth properties		
Citation	WEHI-3 Cytion 400381	
Biosafety level	2	
NCBI_TaxID	10090	
CellosaurusAccession	CVCL_3622	
Receptors expressed	Fc	C3

WEHI-3 | 400381

Viruses

Products	G-CSF	-3	3 IL-3
-----------------	-------	----	--------

Culture Medium	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO ₃	Cytion	820700a
-----------------------	--------------------	------------------------------	--------	---------

Supplements 10% FBS

Subculturing	⁵ / ⁵ 1×10 ⁶ / 3×10	5×10
---------------------	--	------

Fluid renewal 2 3

Freeze medium	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

Thawing and Culturing Cells

- 1.
2. -150°C 3
3. 37°C 40-60
4. 70%
5. 8 15
6. 300 x g 3
7. 10 T25 T25
- 8.

Incubation Atmosphere	37°C, 5% CO ₂
------------------------------	--------------------------

Flask Coating

WEHI-3 | 400381

**Freezing
Procedure**

-78 °C

**Shipping
Conditions**

-78 °C

**Storage
Conditions**

-150 -196 -80 °C

/ **/HLA**

Sterility

PCR