

**HNO41 | 300126**

<b>Description</b>	HNO41	HNSCC	3q23-qter 5p 7p 7q21-q22 8q22.2-qter 9q22-qter 11q24
	HNO41		VEGF PDGF

**Organism**

**Tissue**

<b>Disease</b>	HNSCC
----------------	-------

<b>Age</b>	52
------------	----

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

<b>Citation</b>	HNO41 Cytion 300126
-----------------	---------------------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>CellosaurusAccession</b>	CVCL_D224
-----------------------------	-----------

**HNO41 | 300126**

<b>Culture Medium</b>	DMEM w 4.5 / w 4 L- w 3.7 / NaHCO3 w 1.0	Cytion 820300a
-----------------------	--	----------------

<b>Supplements</b>	10% FBS
--------------------	---------

<b>Dissociation Reagent</b>	Accutase
-----------------------------	----------

<b>Subculturing</b>	PBS T25 3-5 PBS T75 5-10	Accutase T25 1-2 T75 2.5
---------------------	--------------------------	--------------------------

<b>Fluid renewal</b>	2 3
----------------------	-----

<b>Freeze medium</b>	FBS +10% DMSO CM-1 Cytion 800100	CM-1 Cytion 800100
----------------------	----------------------------------	--------------------

<b>Thawing and Culturing Cells</b>	1.		
	2.	-150°C	3
	3.	37°C	40-60
	4.	70%	
	5.	8	15
	6.	300 x g 3	
	7.	10	T25 T25
	8.		

<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>
------------------------------	--------------------------

<b>Flask Coating</b>	
----------------------	--

<b>Freezing Procedure</b>	-78 °C
---------------------------	--------

**HNO41 | 300126**

**Shipping  
Conditions**

-78 °C

**Storage  
Conditions**

-150 -196 -80 °C

/ **/HLA**

**Sterility**

PCR