

**HROG04 | 300932**

**Description** - 2006 CRC

**Organism**

**Tissue** R

**Disease** IV

**Age** 53

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

**Citation** HROG04 Cytion 300932

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_4U39

**Antigen expression** HLA-A02 +, MHC class I -MHC class II -, Beta-microglobulin +, HLA-E +, HLA-G -, MIC A -, MIC-B -, ICAM-1 +, GFAP +, nestin +, vimentin +, S-100 +, GBM +, BTSC +

**Mutational profile** PTENW274L 9p212.3 CDKN2A

**HROG04 | 300932**

<b>Culture Medium</b>	DMEM:Ham's F12 1:1	w: 3.1 g/L	w: 2.5 mM L-	w: 15 mM HEPES	w: 0.5 mM	w: 1.2 g/L NaHCO3	Cytion	8
-----------------------	--------------------	------------	--------------	----------------	-----------	-------------------	--------	---

<b>Supplements</b>	10% FBS
--------------------	---------

<b>Dissociation Reagent</b>	Accutase
-----------------------------	----------

<b>Subculturing</b>	PBS	T25	3-5	PBS	T75	5-10	Accutase	T25	1-2	T75	2.5
---------------------	-----	-----	-----	-----	-----	------	----------	-----	-----	-----	-----

<b>Seeding density</b>	$1 \times 10^4$	/
------------------------	-----------------	---

<b>Fluid renewal</b>	3	5
----------------------	---	---

<b>Freeze medium</b>	50%	+ 40% FBS + 10% DMSO	CM-1	Cytion	800100
----------------------	-----	----------------------	------	--------	--------

<b>Thawing and Culturing Cells</b>	1.					
	2.		-150°C			3
	3.		37°C		40-60	
	4.			70%		
	5.		8	15		
	6.	300 x g	3			
	7.	10			T25	T25
	8.					

<b>Incubation Atmosphere</b>	37°C, 5% CO2
------------------------------	--------------

**HROG04 | 300932**

---

**Flask Coating**

**Freezing Procedure**

-78 °C

**Shipping Conditions**

-78 °C

**Storage Conditions**

-150 -196 -80 °C

/ **/HLA**

**Sterility**

PCR

**HLA**

**A\***: '01:01:01, '05:35  
**B\***: '08:01:01, '51:01:01  
**C\***: '07:01:01, '15:02:01  
**DRB1\***: '03:01:01, '11:01:01  
**DQA1\***: '05:01:01, '05:05:01  
**DQB1\***: '02:01:01, '03:01:01  
**DPB1\***: '01:01:01, '04:01:01  
**E**: '01:01:01