

Product sheet

HROC39 T0 M2 | 300821

Description	-	2006	CRC
Organism			
Tissue	UICC lib	CRC	TNM T4N0M0R0L0V1 G3 Lk(n) + 0 Σ Lk(n) 34
Disease			
Synonyms	HROC39x		
Age	69		
Gender			
Ethnicity			
Morphology			
Growth properties			
Citation	HROC39 T0 M2 Cytion	300821	
Biosafety level	1		
NCBI_TaxID	9606		
CellosaurusAccession	CVCL_1U82		
Protein expression	PTEN		

HROC39 T0 M2 | 300821

Antigen expression Her2/neu + +

Tumorigenic

Viruses SV40 JC/BK HBV HCV HIV

Ploidy status

MSI-status MSS

Mutational profile APCmut, K-Raswt, N-Ras wt, H-Raswt SNP rs12628 at codon 27, B-Rafwt, p53 wt, PIK3CAwt

Culture Medium DMEM:Ham's F12 1:1 w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO₃ Cytion 8

Supplements 10% FBS

Dissociation Reagent Accutase

Doubling time 38

Subculturing PBS T25 3-5 PBS T75 5-10 Accutase T25 1-2 T75 2.5

Seeding density 2×10^4 /

Fluid renewal 3 5

Post-Thaw Recovery

Freeze medium FBS +10% DMSO CM-1 Cytion 800100 CM-1 Cytion 800100

Product sheet

HROC39 T0 M2 | 300821

Thawing and
Culturing Cells

- | | | | | | |
|----|---------|--------|-----|-------|-----|
| 1. | | | | | |
| 2. | | -150°C | | | 3 |
| 3. | | 37°C | | 40-60 | |
| 4. | | | 70% | | |
| 5. | | 8 | 15 | | |
| 6. | 300 x g | 3 | | | |
| 7. | 10 | | | T25 | T25 |
| 8. | | | | | |

Incubation
Atmosphere

37°C, 5% CO₂

Flask Coating

Freezing
Procedure

-78 °C

Shipping
Conditions

-78 °C

Storage
Conditions

-150 -196 -80 °C

/ /HLA

Sterility

PCR