

OKT 3 | 300499

**Description** P3x63Ag8.U1 Muromonab-CD3

**Organism**

**Tissue**

**Applications**

**Synonyms** OKT3

**Breed/Subspecies** CAF1 x BALB/c

**Morphology**

**Cell type** B

**Growth properties** /

**Citation** OKT 3 Cytion 300499

**Biosafety level** 1

**NCBI\_TaxID** 10090

**CellSaurusAccession** CVCL\_2665

**Products** T T CD3

OKT 3 | 300499

<b>Culture Medium</b>	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3 Cytion	820700a
-----------------------	--------------------	-------------------------	---------

<b>Supplements</b>	10%	FBS
--------------------	-----	-----

<b>Subculturing</b>	<sup>5</sup> /	1×10
---------------------	----------------	------

<b>Freeze medium</b>	FBS +10% DMSO	CM-1 Cytion	800100	CM-1 Cytion	800100
----------------------	---------------	-------------	--------	-------------	--------

<b>Thawing and Culturing Cells</b>	1.			
	2.	-150°C		3
	3.	37°C	40-60	
	4.		70%	
	5.	8	15	
	6.	300 x g	3	
	7.	10	T25	T25
	8.			

<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>
------------------------------	--------------------------

<b>Flask Coating</b>
----------------------

<b>Freezing Procedure</b>	-78 °C
---------------------------	--------

<b>Shipping Conditions</b>	-78 °C
----------------------------	--------

