

Product sheet

5A w 3.0 / w w 2.0 w 2
.2 / NaHCO3 | 820200a

5A
G LGR5 HCT116 RKO FET CBS HCT116b TENN
5A
5A Eagle
5A 5A

5A 3.0 g/L 2.0 mM 2.2 g/L

- pH =7.2±0.02 20-25°C
-
- +2°C +8°C +37°C
- 37°C
-
- 6 8

		/
2		132.00
		97.67
		400.00
		6,460.00
		504.00
D(+)-		3,000.00
		0.50
		600.00
		11.00

5A w 3.0 / w w 2.0 w 2
 .2 / NaHCO₃ | 820200a

L-				13.36
L-	x			42.14
L-	x			45.03
L-				19.97
L-	x	x		31.75
L-				219.15
L-				22.07
				7.51
L-	x	x		20.96
L-				19.67
L-				39.36
L-				39.36
L-	x			36.54
L-				14.92
L-				16.52
L-				17.27
L-				26.28
L-				17.87
L-				3.06
L-				26.10

Product sheet

5A w 3.0 / w w 2.0 w 2
.2 / NaHCO3 | 820200a

L-	17.57
	1.00
	0.56
D(+)-	0.20
D-	0.20
	5.00
	10.00
	36.00
	0.50
	0.50
	0.50
	0.50
	0.20
	0.20
B12	2.00