

L-WRN | 300641

Description	L-WRN	L	L	L-WRN	Wnt3a R-3	Noggin
		L-WRN	Wnt	Wnt	L-WRN	Wnt3a R-3
	L-WRN				Wnt	

Organism	
Tissue	

Applications	
---------------------	--

Breed/Subspecies	C3H/An
-------------------------	--------

Age	100
------------	-----

Gender	
---------------	--

Morphology	
-------------------	--

Growth properties	
--------------------------	--

Citation	L-WRN Cytion 300641
-----------------	---------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_DA06
-----------------------------	-----------

GMO Status	GMO-S1 NIH-3T3 L-WRN Wnt3a R-spondin-3 Noggin SV40 DNA hph Tn5-neo
-------------------	--

L-WRN | 300641

Freezing Procedure -78 °C

Shipping Conditions -78 °C

Storage Conditions -150 -196 -80 °C

/ /HLA

Sterility PCR