

NCI-H1395 | 305118

<b>Description</b>	NCI-H1395	(NSCLC)	55	15	1986
--------------------	-----------	---------	----	----	------

**Organism**

**Tissue**

**Disease**

<b>Synonyms</b>	NCI-H1395 H-1395 NCIH1395
-----------------	---------------------------

<b>Age</b>	55
------------	----

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

<b>Citation</b>	NCI-H1395 Cytion 305118
-----------------	-------------------------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>CellosaurusAccession</b>	CVCL_1467
-----------------------------	-----------

NCI-H1395 | 305118

<b>Culture Medium</b>	RPMI 1640 w 2.0 mM	w 2.0 g/L NaHCO3	Cytion	820700a
-----------------------	--------------------	------------------	--------	---------

<b>Supplements</b>	10% FBS
--------------------	---------

<b>Dissociation Reagent</b>	Accutase
-----------------------------	----------

<b>Subculturing</b>	PBS	T25	3-5	PBS	T75	5-10	Accutase	T25	1-2	T75	2.5
---------------------	-----	-----	-----	-----	-----	------	----------	-----	-----	-----	-----

<b>Fluid renewal</b>	2	3
----------------------	---	---

<b>Freeze medium</b>	FBS +10% DMSO	CM-1	Cytion	800100	CM-1	Cytion	800100
----------------------	---------------	------	--------	--------	------	--------	--------

<b>Thawing and Culturing Cells</b>	1.										
	2.		-150°C					3			
	3.		37°C			40-60					
	4.			70%							
	5.			8	15						
	6.	300 x g	3								
	7.	10				T25		T25			
	8.										

<b>Incubation Atmosphere</b>	37°C, 5% CO2
------------------------------	--------------

<b>Flask Coating</b>	
----------------------	--

<b>Freezing Procedure</b>	-78 °C
---------------------------	--------

NCI-H1395 | 305118

**Shipping  
Conditions**

-78 °C

**Storage  
Conditions**

-150 -196 -80 °C

/ /HLA

**Sterility**

PCR